



## Patient Preparation /Collection Instructions

- x Cleanse skin or mucosal surfaces. For closed wounds and aspirates, disinfect as for a blood culture collection with 2% chlorhexidine or 70% alcohol followed by an iodine solution. Remove iodine with alcohol prior to specimen collection. For open wounds, debride (if appropriate), and thoroughly cleanse with sterile saline prior to collection. Sample viable infected tissue, rather than superficial debris.
- x Tissue biopsies should be collected from areas within and adjacent to the area of infection. (large enough tissue samples should be collected to perform all of the tests requested)

## Performance

### Days Performed

Daily; Monday - Sunday

### Report Available (TAT) (Once received at MDL)

4 - 6 days

### Specimen Retention Time

7 days

### Method Description

- x Conventional aerobic and anaerobic bacterial culture technique with selective and nonselective media.
- x Identification methods (when appropriate) may include any of the following: conventional biochemical testing, matrix-assisted laser desorption/ionization time-of-flight (MALDI-TOF) mass spectrometry, and commercial identification panels.
- x Susceptibility testing (when appropriate) may include minimal inhibitory concentration (MIC) (broth microdilution or gradient strip diffusion) or disk diffusion.

### Reference Values

No growth.

## Cautions

- x Antibiotics administered prior to sample collection may negatively affect the recovery of organisms associated with infection. Preferably collect specimen prior to initiation of therapy and only from wounds that are clinically infected or deteriorating or that fail to heal over a long period.
- x Many wound infections are polymicrobial and the isolation of an organism in culture may or may not correlate with infection of the wound.